

Near Field Scanning Optical Microscopy (NSOM)

PH 575 Presentation

Ram Prakash (210260042), Aditya Mehta (210260032)

IIT Bombay

October 24, 2024



Outline

- 1 Introduction
- 2 Instrumentation
- 3 Working Principle
- 4 Applications
- 5 Limitations
- 6 Recent Developments
- 7 Conclusion

What is NSOM?

- Near Field Scanning Optical Microscopy (NSOM), also known as Scanning Near-field Optical Microscopy (SNOM), is a microscopy technique that surpasses the diffraction limit of light.
- It uses a subwavelength probe to scan the sample at a very close distance (within the near field).
- Provides high spatial resolution beyond conventional optical microscopes.
- Traditional optical microscopy is limited by diffraction to about half the wavelength of light (≈ 200 nm for visible light).
- NSOM breaks this barrier, achieving resolution in the range of 20-100 nm.
- Enables the study of biological samples, semiconductors, and nanostructures at a higher resolution.

NSOM Setup




FIG. 2. Standard SNOM setup consisting of (a) an illumination unit, (b) collection and redistribution unit, and (c) a detection module.

Figure: Image taken from [J. Chem. Phys. 112, 7761–7774 \(2000\)](#).

1. Optical Unit: Components

- **Laser Source:**

- Provides coherent light necessary for near-field measurements.
- Common types: Semiconductor lasers.

- **Optical Fibers:**

- Transmits laser light to the SNOM probe with minimal loss.
- Can be single-mode or multi-mode, depending on application.

- **SNOM Probe:**

- Equipped with a sharp tip (typically made of metal) for near-field interaction.
- Small aperture size (10-100 nm) to achieve high spatial resolution.

- **Shear Force Sensor**

- Used for feedback mechanism

2. Scanning Stage Unit

- **High-Precision Scanning Stage:**

- Allows for three-dimensional movement of the SNOM probe.

- **Feedback Mechanisms:**

- Ensures consistent distance from the sample for accurate measurements.
- Moves the probe relative to the sample in a controlled manner.
- Maintains optimal height for accurate near-field measurements using real-time feedback.
- Capable of rapid scanning over large areas while maintaining high precision.
- Ensures the probe tip remains within the near-field zone for effective signal collection.

Feedback Mechanism in the Scanning Stage

- **Purpose:**

- Maintain optimal probe-sample distance for accurate measurements.

- **Types of Feedback Mechanisms:**

- **Tunneling Current Feedback:**

- Monitors tunnelling current to keep a constant distance.

- **Shear Force Feedback:**

- The fiber probe is vibrated at its mechanical resonance parallel to the sample surface
- Amplitude and phase of the fiber's oscillation are monitored using a suitable displacement sensor
- As the probe approaches the sample (0–20 nm), shear forces cause a detuning of the resonance frequency, resulting in a decrease in amplitude and a phase shift

- **Benefits:**

- Ensures stability and enhances resolution.
- Improves accuracy by compensating for position deviations.

3. Detection and Data Acquisition Unit

- **Components:**

- Detectors (e.g., PMTs)
- Signal processing electronics
- Data acquisition systems

- **Function:**

- Collects light emitted/scattered from the sample
- Converts intensity into electrical signals for analysis and image formation

- **Transmission Mode:**

- The light is transmitted through the sample
- Useful for thin samples or transparent materials
- Provides information about the internal structure and properties
- Typically used in studies of layered materials and biological samples

- **Reflection Mode:**

- Detects light reflected from the sample surface
- Ideal for opaque materials where light cannot pass through

Detection Modes in SNOM (Cont.)

- **Reflection Mode:**

- Provides surface topography and optical property information
- Can reveal features like roughness and composition differences

- **Collection Mode:**

- The probe collects scattered light from the sample surface
- Sensitive to local features and can detect changes in optical properties
- Useful for studying localized phenomena such as fluorescence and surface plasmon resonance
- Enables high-resolution imaging of nanostructures and biomolecules

- **Comparison of Modes:**



- Each mode has specific applications and advantages depending on the sample type
- Modes can be combined for multi-modal imaging, enhancing the depth of analysis

A thought Experiment

Suppose I have a glass prism, oriented as shown, with a laser beam undergoing total internal reflection from the internal surface.

Question: at what point do I see a transmitted beam?

Answer: when the prisms are close together, but not yet touching!



Evanescent Waves and Optical Resolution

- **Evanescent Waves:**

- Generated when light undergoes total internal reflection (TIR) at an interface.
- Characterized by an exponentially decaying amplitude in the near-field region.
- Electric and magnetic fields in a harmonic wave can be described as:

$$E = E_0 e^{i(k \cdot r - \omega t)}, \quad H = H_0 e^{i(k \cdot r - \omega t)}$$

- **Total Internal Reflection (TIR):** Occurs when light incident on an interface exceeds the critical angle:

$$\sin \theta_c = \frac{n_2}{n_1}$$

where n_1 and n_2 are the refractive indices of the media, and θ_c is the critical angle.

Evanescent Waves and Optical Resolution

- **Evanescent Wave Equation:** For angles greater than the critical angle, the transmitted wave becomes evanescent, with an exponentially decaying amplitude:

$$E_t = E_0 e^{-\alpha y}$$

where $\alpha = k_t \sqrt{\sin^2 \theta - \left(\frac{n_2}{n_1}\right)^2}$ is the decay constant.

- **Skin Depth:** The skin depth δ in the visible spectrum (e.g., for metals) is given by:

$$\delta = \frac{1}{\alpha}$$

Typical values for the skin depth in good conductors are around 10 nm in the visible spectrum.

Working Principle of NSOM

- NSOM captures near-field evanescent waves, which decay quickly over distance but carry high spatial frequency information.
- A nanoscale probe is placed extremely close to the sample, within a few nanometers, to detect these near-field interactions.
- The probe interacts with the sample, converting damped waves into propagating waves, which are then detected to form a sub-wavelength resolution image.
- The probe in NSOM is often a tapered optical fiber coated with metal, with apertures as small as 10-100 nm. For aperture SNOM, the resolution is determined by the probe's aperture, while apertureless SNOM relies on scattering at a sharp tip.
- In a nutshell, Light interacts in the near-field region, allowing for the detection of light scattered, absorbed, or phase-shifted due to sample interactions, with high-resolution images formed from these localized fields

Probe and Light Interaction




Figure: Probe and Light interaction

- High-resolution imaging of nanostructures and biological samples, with applications in detecting defects, dopants, and molecular structures in semiconductors.
- Combining SNOM with Raman spectroscopy enables high spatial resolution (down to 30 nm) for studying molecular structures and phases.
- Single molecule detection is easily achievable.
- Dynamic properties can also be studied easily at sub-wavelength scales.

Applications of NSOM

SNOM is presented as a powerful tool for the subdiffraction limited chemical and structural characterization of plant cell walls. It provides a detailed knowledge on the macromolecular assembly of plant cell walls.




Figure: Analyzing Plant Cell Wall Ultrastructure

Limitations of NSOM

- Probe fabrication, particularly achieving a small and smooth aperture (10-100 nm), is difficult and affects resolution.
- Sample-probe interaction can lead to image artifacts, especially if the tip-sample distance isn't precisely controlled.
- Slow scan rates are a limitation, and the requirement for highly skilled operators further hinders widespread adoption.
- Not conducive for studying soft materials, especially under shear force mode

Recent Developments in NSOM

- Development of aperture-less NSOM for higher resolution.
- Use of plasmonic effects to enhance near-field interactions.
- Hybrid NSOM systems combining AFM for enhanced resolution and contrast.
- Integration of NSOM with spectroscopic techniques for chemical analysis.

- NSOM has proven to be an indispensable tool for studying nano-scale structures that are beyond the resolution of traditional optical microscopy. By capturing near-field light interactions, NSOM offers unmatched spatial resolution for a wide range of applications in material science, biology, and nanotechnology.
- While challenges such as probe fabrication and scan speeds remain, ongoing advancements in probe design, feedback mechanisms, and detection methods are continually improving NSOM's accuracy and resolution

- Hecht, B. et al. (2000). Scanning near-field optical microscopy with aperture probes: Fundamentals and applications. *Journal of Chemical Physics*, 112(18), 7761–7774.
- Sobat, Z., Hassani, S. S. (2014). An overview of scanning near-field optical microscopy in characterization of nano-materials. *International Journal of Nano Dimension*, 5(3), 203-212.
- Ajay Pratap, (2020). Scanning Near Field Optical Microscopy. *Module VI: Nano Material*.